

Table S1 – RNA extraction protocol

All materials are to be treated RNasefree.

- (1) Lyse cells with Trizol (1 ml per 1×10^6 cells) and sonicate for 5 to 10 min.
- (2) Add 0.2 ml chloroform/ml.
- (3) Shake vigorously for 30 s and incubate 5 min at RT.
- (4) Centrifuge 15 min at 12000 x g, 4 °C.
- (5) Transfer aqueous phase into new tube.
- (6) Add 0.5 ml isopropyl alcohol, 0.12 ml 5M ammonium acetate, 10 μ l of 0.5 mg/ μ l glycogen. Mix.
- (7) Incubate 20 min at -20 °C.
- (8) Centrifuge 30 min at 12000 x g, 4 °C.
- (9) Wash pellet with 75% EtOH and resuspend in 85 μ l H₂O.
- (10) Perform DNase digest.
- (11) Perform phenol/chloroform/isoamyl alcohol (25:24:1) extraction.
- (12) Add 0.1 v/v 3M sodium acetate (pH 5.5) and 2.0 v/v EtOH.
- (13) Incubate at -20 °C overnight.
- (14) Centrifuge 30 min at 12000 x g, 4 °C.
- (15) Resuspend in appropriate amount of H₂O.