

Intonation Analysis System

- Users Guide -

Stand: 31.10.1981

Verf.: Peter Wittenburg, Gunter Senft

Inhalt:

1. intention
2. global layout
3. hardware components/cabling
4. software components/parameters
5. procedure / example

1. Intention

The intonation analysis system serves to calculate resp. display the suprasegmental parameters of a given piece of natural speech.


These are: - basic frequency (f_0) \Rightarrow pitch
- amplitude \Rightarrow intensity
- durations \Rightarrow time pattern

It also can serve to test pitch detection algorithms based on short term Fourier representations and to display series of these Fourier representations (f-domain representations) to detect other parameters.

The set-up was chosen so that nearly real time performance is achieved, i.e. that the user will get an output without any large delay. So the system works rather efficiently both in "interactive" sessions with short speech material and with large amounts of speech to be analyzed.

2. Global layout


a) Hardware



The speech can be entered via Micro or TR. It passes some fast, special purpose electronics, so that a conversion from the original speech into series of spectral representations and a coding of the intensity. Also synchronization purposes pulses can be entered into the information, which can be seen later on on the graphic output.

On the computer, three programs serve for data acquisition, analysis and output. The output is first generated on a graphic display and then a copy can be made on paper.

b) Software



The software is set up by starting the indirect command file DL:INT , which resides on the disc DL:Intonation.

The acquisition is done by starting the indirect command file INTIST . This program creates a file `.SPC` which contains series of spectra and a file, `.INT` which contains the sampled intensity values. The derivation of the basic frequency is performed by the program FZE. FZE checks first if the data is acquired correctly and then calculates per spectrum the F_0 -value which is stored sequentially in a file `.FZE`. Therefore some analysis parameters have to be defined. The program OUT performs the output of both F_0 and intensity (if given also sync marks are drawn). Some parameters defining the form of the output have to be entered.

Figure 1 shows the correct cabling of the system. The Intonation Analysis System (IAS) can only work properly, if the cabling is done correctly. This is, of course, a rather trivial statement, but experience has shown that troubles start when people try to do the cabling. Therefore make sure that:

- 1st the envelope generator (2) is connected with
 - the attenuator (3)and with the computer via the
 - AD 11K cable
- 2nd the attenuator (3) is connected with the preamplifier (4)
(and, of course, with the envelope generator (2))

3rd the preamplifier (4) is connected with the tape-recorder and with the spectrum analyzer (and, of course, with the attenuator (3)).

4th the system is connected with the computer:

therefore check at the computer whether the two cables (label 125) are connected with the DR 11 B via the front side connectors.

If you have to do that connection at the computer, handle the cables very carefully!

5th the computer has to be loaded with the INTONATION DISC

If you have to load the computer with the intonation disc


- ask someone from the technical group, how and where you can load the disc.
- press the load button
- wait till load button lights up
- change the discs, that is put in your intonation disc and make sure that it fits correctly in the machine.
- press the load button
- wait till the "ready" button lights up.

Figure 1

2. hardware/cables

check carefully

- 1 = Pitch Module
- 2 = Envelope generator
- 3 = Attenuator
- 4 = Pre-Amplifier



The whole setup of the hardware should normally be chosen as indicated on the photo (see photo 1).

If you have done the cabling, switch on all the units needed for the intonation analysis; these are:

- the tape recorder
- the spectrum analyzer
- the scope
- the terminal
- and
- the hard copy unit

Please make sure after having finished your analysis that all these systems are switched off before you leave the room.

If all the units remain switched on for a very long time (e.g. over night) they may cause severe damage to the units.

The most important point is to adjust all auditive levels so that the full dynamic range of all components is optimally used, i.e. test for a sentence that the amplifiers resp. attenuators are adjusted so that the overflow/overload lamps on spectrum analyzer, amplifier, and envelope generator just do not come up.

Adjusting the amplifier needs turning up or down the volume control knob. Adjusting the attenuator needs setting up or down the number indicated at the top of this unit.

A correct adjustment of all auditive levels is the basis of a good intonation analysis. To do it properly needs some experience with the IAS. So do some tests, "play" with the system, and do not give up too soon.


On the amplifier (4) you can chose whether you want to use the microphone or the tape recorder as "speech source", i.e. check the position of the switch before you start.

On the spectrum analyzer you should check whether the settings are exactly as indicated on the photo, except the input attenuation which has to be adjusted by the user.

When switching the spectrum analyzer on you have to position the cursor until you have got a number between 500 and 1000 with the last digit being a zero (examples 500, 680, 970, ...; not allowed: 791).


This number divided by 10 is used in the software as the parameter N!!

The normal value should be 500, so $N = 50$.



The oscilloscope gives you a kind of first hint about your input by its visual display.

However, do not take this display too serious. You can adjust the base line of the display by turning the two variable knobs on the right of the unit.



4. Software/parameters

Software operation includes 7 steps, which have to be done by the user:

- login and mount disc
- run command file $\text{\textcircled{a}}$ INTONAT
- run acquisition program $\text{\textcircled{a}}$ INTST
- run analysis program FZE
- run output program OUT
- delete files (if you do not need them anymore)
- run command file $\text{\textcircled{a}}$ STOPINT

The software exists out of 5 components:

1. a start-up command file INTONAT.CMD
2. a data acquisition indirect command file $\text{\textcircled{a}}$ INTST
3. a basic frequency detection program FZE
4. an output program OUT
5. a shut -down command file STOPINT.CMD

The whole software resides on DL:INTONAT

1. INTONAT.CMD

Drivers are loaded and tasks are installed.

Program INT is started.

2. $\text{\textcircled{a}}$ INTST

INT samples data from spectrum analyzer and envelope generator. Spectral data is stored in a file .SPC while the intensity data (sampling frequency 100 Hz) is stored in a file .INT.

INT at first asks for a file specification: {DL:INT } < name >

use DL: if you do not have to much sentences to analyze at one step

use MT: if you have several sentences to be analyzed at one step

Name can be chosen by the user (not more than 7 letters).

then INT asks for three parameters: NSEC, N, RANGE (all integers)

NSEC you have to input an estimate of the length of your stimuli (in sec)

N you have to divide the number indicated on the spectrum analyzer by ten.

RANGE band setting on spectrum analyzer (normally 5 letters)

INT then opens a file, if opening was successful INT asks whether everything is OK, where you can answer Y or N (instead of Y also return) INT then prepares everything and prompts with "RETURN TO START SAMPLING". It starts immediately sampling after having typed in "return". If you work with the tape-recorder, type first return, then start the tape.

3. FZE

FZE checks the correctness of the data stream (synchronity check), decodes the bit pattern and calculates for every spectrum the basic frequency F.

It at first asks for a file specification {DL:/MT} < name >
So you just enter the same as you have entered in the program INT

FZE then opens the files: {DL:/MT:} < name > .SPC } old
DL: < name > .INT
DL: < name > .FER new

From file .INT it reads from the first block the parameters N, RANGE as default values.

FZE then asks whether you want to input parameters; if you type N, it uses the default values N = 50, RANGE = 5, PAUSE = 1000, FTHR = 8, NTHR = 20, YTHR = 1.7, XTHR = 5.0.

These settings turned out to work pretty well for a wide range of material.

Meanings of some parameters:

.FTHR serves as an ideal high pass filter, since FTHR = 8 means that all frequency components smaller than 80 Hz are ignored. For female child speech you can make FTHR even higher.

.NTHR serves as a noise rejection level (in voltage), so depending on the sound level you can vary this parameter.

.YTHR and XTHR are parameters used by the pitch detector 1.7 and 4.0-6.0 turn out to deliver good results.

.PAUSE serves to identify pauses, i.e. when the number of spectra with all values < NTHR exceeds PAUSE than a pause is indicated and calculation starts again when a value > NTHR.

4. OUT

OUT is the program which performs the output of intensity, F_0 and sync pulses (if present). With help of parameters you can vary the format of the output.

OUT at first asks for a file name only: < name > since the device always is DL:. It opens two files: DL: < name > .INT }old

DL: < name > .FZR

It then asks if you want to have the intensity curve in dB (D) or metric scale (V). If you don't want to use the standard default values, you have to enter new parameter-values:

- POINTS/STEP is a factor which just stretches the output format on the display.

.DELTAT determines the time resolution of your output, i.e. if DELTAT = 10 then you get for each 10 ms a value of both I and F_0 (if necessary a simple interpolation is done).

.FTHR and PAUSETHR again serve to identify pauses to reduce the amount of output

.FLOW and FHIGH determine the F_0 sealing.

.SCALE-STEP determine the sealing on the time axis in msec.

5. Procedure/Examples

- Turn on power of all units at least 5 minutes before you start working!!
i.e. switch at the side of the rack system, switch of hardcopy unit,
switch of spectrum analyzer and scope.
- press erase on display
- check all settings, as indicated in this documentation
- check all cabling, as indicated in this documentation
(if the lights are still on on the Pitch Module, the cable of this
computer is not connected properly)
- check the microphone (TR switch)
- adjust the audio levels so that no overflow occurs
- check, if you want to use a Mag Tape, if it is installed
- then login, mount the disc and start the software (INTONAT)

operation
- when finished stop the software (STOPINT) (BYE)
- switch power off on all units

type
HEL 7,1/PETER

RSX-11M BL26 MULTI-USER SYSTEM

GOOD AFTERNOON
15-NOV-81 15:35 LOGGED ON TERMINAL TT4:

Welcome to RSX-11M V3.2 timesharing

* ***** You are logged onto the PDP 11/55, our mainframe!!! *****

>MOUNT ** VOLUME INFORMATION **

CLASS = FILES 11

DEVICE = DL0:

LABEL = INTONATION

UIC = [1,1]

VOL PRO = [RWED,RWED,RWED,RWED]

FILE PRO = [RWED,RWED,RWED,R]

CHAR = []

ACP NAME = F11ACP

>

type

———— check if Label "Intonation" is printed

type

```
@INTONAT
>UNL SP:
UNL -- DEVICE SP0: NOT IN SYSTEM
>UNL XB:
>UNL MT:
>LOA LA:
LOA -- DRIVER ALREADY RESIDENT
>LOA XB:
>LOA MT:
>REM INT
>REM FZE
>REM OUT
>INS DL:INT/TASK=...INT/PRI=200./CKP=NO
>FIX INT
>INS DL:FZERA/TASK=...FZE
>INS DL:INTOUT/TASK=...OUT
>@ <EOF>
>
```

type

@INTST
>REM INT
>INS INT/TASK=...INT/PRI=200.
>INT

type

SPECIFY DEVICE AND FILENAME (NO EXTENSION) : DL:SATZ1
INTFIL : DL:SATZ1.INT } two files are created
FOFIL : DL:SATZ1.SPC

type

SPECIFY NSEC, N AND RANGE : 4,50,5

type

DMA PARAMETERS :
NBUF : 9 NSEC : 4 BUFSIZ : 4096

type

LPA PARAMETERS :
NBUF : 1 NSEC : 10 BUFSIZ : 1024

type

ALL PARAMETERS OK ?

type

Return
RETURN TO START SAMPLING

type

Return

(or DL: Refer 1 etc)

disc name to be chosen by the user

4

nr. of seconds

50

(500/10)
nr. of points

5 (kHz)

setting of range

* Return means: type "Return" key on the terminal

DMADON
 DMBFDN
 DMADON
 DMBFDN
 DMADON
 DMBFDN
 DMADON
 DMBFDN
 DMADON
 DMBFDN
 DMADON
 DMBFDN
 DMADON
 DMBFDN
 DMADON
 DMBFDN
 DMADON
 DMBFDN
 DMADON
 DMBFDN
 DMADON
 DMBFDN
 DMADON
 DMBFDN
 DMADON
 DMBFDN
 DMADON
 DMBFDN
 DMADON
 DMBFDN
 DMADON
 ALL DMA BUFFERS FILLED
 DMBFDN
 LPADONE
 LP-DK-DONE
 ALL REQUESTS WELL DONE
 ...INT -- STOP END OF INTONATION SAMPLING
 > < EOF >
 >

input in progress

end of input ; stop speaking

input was successful

type
type

FZE
SPECIFY DEVICE AND FILENAME (NO EXTENSION) : DL:SATZ1

make sure that you use the same name as in

DL:SATZ1.INT
DL:SATZ1.FZR new file is created

DL:SATZ1.SPC
PARAMETER INPUT WISHED ? Y

⊙ INTST

N,RANGE,PAUSE : 50,5,1000

nr. of points in spectrum, range setting, nr. of points after which a pause is identified

FTHR,NTHR,YTHR,XTHR : 8,20,1.7,5.0

START AT BLOCKNUMBER : 17

NEW BUFFER READ
NEW BUFFER READ
NEW BUFFER READ
NEW BUFFER READ
NEW BUFFER READ
NEW BUFFER READ
NEW BUFFER READ
NEW BUFFER READ

8 : (= 80 Hz)
lowest expected frequency

-20 : metric noise rejection level
use a higher level, if your tape is noisy

1.7 } algorithm dependent variables
5.0 }

EOF REACHED
FZER TERMINATED *****

type
type
type
type
type
type
type
type

```

OUT
INPUT FILENAME : SATZ1
INTENSITY FILE : DL0:SATZ1.INT
FZERO FILE : DL0:SATZ1.FZR
INT.FILE IN VOLT- OR DB-SCALE(U/D) : D
PARAMETER CHANGE WISHED(Y/N) : Y
INPUT POINTS/STEP, DELTAT : 4, 10
INPUT FTHR AND PAUSETHR : 1, 1000
INPUT FLOW, FHIGH, SCALE-STEP : 50, 250, 100
AUTOMATIC COPY(Y/N) : N
MEAN(Y/N) : Y
  
```

name

D = logarithmische Output
 V = metrischer Output

if N; default values are taken
 4; stretches (use then 5, 6, or any higher integer) or compresses (use 3 or any lower integer) the output on paper

10: 10 ms per drawn point f₀

if Y automatic hardcopy is made if you have a number of sentences, pauses are identified and not drawn on paper

mean value is calculated and drawn as dotted line

50-250 = upper or lower boundary for frequency scaling

high voices (female voices) most often need a higher boundary than 250


100 = 100 ms scaling on time axis

the terminal then shows the analysis in such a display:


FILE : PETER1 DELTAT : 10 ^{ms}
DELTA F : 6.23 ^{ms} DELTA I : 10.00 ^{ms}
SAMPL.FREQ. : 100 SCALE : 100 ± 100ms
PAGE : 1

page indication

envelope resp. intensity is sampled with 200Hz



P e t e r e n z a t e n i n h e t h e k




mean value of F0

TYPE


Type Y if you want a hardcopy

FILE : GUNTER1 DELTAT : 10
DELTA : 6.23 DELTAI : 10.00
SAMPL.FREQ. : 100 SCALE : 100
PAGE : 1



nanche leste vede tief

FILE : GUNTER1 DELTAT : 10
DELTA F : 6.23 DELTA I : 10.00
SAMPL.FREQ. : 100 SCALE : 100
PAGE : 2



... lobe Stimm.e barbe,

if the data are successfully sampled in INT you can always go back either to FZE or to OUT and change your parameters there to improve your output

END OF FILE, NEW FILE WISHED(Y/N) : N
...OUT -- STOP END OF OUTINT

- > if you see content with your output,
- > delete the created file by typing:

> PIP SATZ1.*:*/DE

> blank!

> @INTST

Return (before typing "return" key, check your deletion command!)

if you want to go on with your analysis and proceed as indicated above - (id est: DL: Satz 2 etc)

type

type

type

END OF FILE, NEW FILE WISHED(Y/N) : N
...OUT -- STOP END OF OUTINT

>PIP SATZ10.*;*/DE

>

>

>

>BYE

if you have finished your work

type

HAVE A GOOD AFTERNOON

15-NOV-81 16:13 TT4: LOGGED OFF

>

DMO -- TT4: DISMOUNTED FROM DL0: *** FINAL DISMOUNT ***

SWITCH OFF ALL UNITS !
=====

If you have any problems please contact Peter (Tel: Institute: ext. 122

Private: 226066)

We are also looking forward to any suggestions for improving our system. Check, for instance, whether you think it useful to get with your printed output some information like:

- \bar{X} of F_0
- Ranges of F_0
- s d
- or on filters and how they were used to suppress noise.

Have fun and success with your intonation analyses.

Gunter and Peter