

Cycloparaphenylenes and their Catenanes – Complex Macrocycles Unveiled by Ion Mobility Mass Spectrometry

Wen Zhang, Ali Abdulkarim, Florian E. Golling, Hans Joachim Räder*, and Klaus Müllen*

Abstract: We analyze the insoluble product mixture obtained from cycloparaphenylene (CPP) synthesis via Suzuki coupling and reductive aromatization. Traditional mass spectrometry suggests a homologous series of macrocycles with 12 to 84 phenylene units. Ion mobility mass spectrometry (IMMS), however, unravels an unexpected complexity of isomers with identical chemical formula, but different topologies. Whereas macrocycles containing up to 30 phenylene units show only one structure, the homologue with 36 phenylene units forms at least four different isomers with significant molecular size differences. They can be assigned to catenanes composed of CPPs with 2x18 and 12+24 phenylene units together with the ordinary [36]CPP macrocycle. Most likely, a trefoil knot of the CPP with 36 moieties is also present. For the first time, catenanes can be elucidated in a simple reaction mixture by analyzing their ions in the gas phase which lies beyond the scope of traditional analytical methods.

Cycloparaphenylenes (CPPs) are attractive model compounds as they constitute the smallest elements of carbon nanotubes (CNTs) and “seed particles” for a bottom-up fabrication of CNTs.^[1] Furthermore, their electronic and optoelectronic properties qualify CPPs as active components of molecular electronics and photovoltaics.^[2] CPPs with sizes from 5 to 16 and 18 phenylene rings have been synthesized.^[3] With increasing number of phenylenes, different CPPs are hardly separable *via* conventional methods and have escaped further analysis.

Similar to the first pioneering synthesis of [n]CPPs, our synthetic approach (Scheme 1) is based on Suzuki coupling of 1,4-syn-dimethoxy-2,5-cyclohexadiene derivatives **1** and **2** as “kinked” precursors.^[4] The molecules are prone of macrocyclization as shown in previous works.^[5] We extended the reaction time to 96 hours to ensure complete conversion and introduced phenyl groups as end-cappers to achieve a better distinction of cyclic and linear products. The resulting mixture **3** was separated from insoluble material which we attribute to the formation of acyclic oligomers. Reductive aromatization of the cyclic products **3** by treatment with sodium naphthalenide afforded a mixture of fully conjugated CPPs. Their insolubility prevents separation by chromatographic methods and characterization by solution NMR or single crystal X-ray diffraction.

Herein, we elucidate the complex mixture using matrix-assisted laser desorption ionization (MALDI) together with ion mobility mass spectrometry (IMMS). These methods combine three features: i) solvent-free sample preparation of insoluble analytes, ii) molecular size dependent separation by IMS, and iii) tandem mass spectrometry. The analysis points toward homo- and heterocatenanes where only “macro-mono-cycles” would have been expected.

High resolution MALDI MS using our solvent-free sample preparation^[6] (Figure 1) discloses signals up to 6390 Da. The mass difference between neighboring monoisotopic peaks is 456.19 Da, corresponding to the molecular weight of six *para*-phenylene groups. The detected absolute masses correspond exactly to multiples of 456.19 Da, indicating the absence of end groups. Thus, the oligomer peaks can be assigned exclusively to CPPs consisting of equal numbers of the initial A-A and B-B type building blocks with ring sizes from 12, 18, 24 ... up to 84 phenylene rings. A CPP with 6 phenylene units is not detected which can be rationalized by its high ring strain. The absence of odd numbered repeat units, such as oligomers with 3, 9, 15, etc. phenylene rings, further supports the assignment to cyclic structures, because those are not possible from the underlying chemistry. The homologous series of molecules is denoted as P_n where “n” represents the number of phenylene rings. High resolution and exact mass measurements enable the assignment of elemental compositions for each oligomer. As shown for the signal at 2737.1323 Da (insert in Figure 1), the chemical formula is C₂₁₆H₁₄₄ and the experimental isotopic distribution fits well to theory, thus matching a CPP consisting of 36 phenylene units, here denoted as P₃₆.

IMS can separate gas phase ions of equal molecular weight due to their molecular shape differences, represented by their collision cross-section (CCS) values. Even isomeric species can be distinguished, thereby overcoming the intrinsic limitation of mass spectrometry to separate ions exclusively due to their





Figure 1. MALDI mass spectrum of the insoluble CPP product mixture.

[a] Dr. W. Zhang, A. Abdulkarim, Dr. F. E. Golling[†], Dr. H. J. Räder, Prof. Dr. K. Müllen
Max Planck Institute for Polymer Research
Ackermannweg 10, 55128, Mainz (Germany)
E-mail: raeder@mpip-mainz.mpg.de
muellen@mpip-mainz.mpg.de

[†]Present address: Covestro AG, 51365, Leverkusen, Germany

Supporting information for this article is given via a link at the end of the document.

Scheme 1. Synthetic route towards cycloparaphenylenes (CPPs).



mass-to-charge (m/z) ratio.^[7] So far, IMS has been used to reveal the molecular architectures of a series of mechanically interlocked polyrotaxanes^[8], catenanes^[9] and polymers^[10].

Figure 2 shows the ion mobility spectra of oligomers up to P_{36} . The measurements are carried out by using 10V collision energy prior to ion mobility separation in order to dissociate non-covalently bonded aggregates of small CPPs which also contribute to a small extend to the signals P_{24} , P_{30} and P_{36} (Figure S1). Therefore, the ion mobility traces in Figure 2 are exclusively due to covalently and mechanically bonded isomers. The ion mobility spectra of the small oligomers P_{12} , P_{18} , P_{24} and P_{30} show only one peak in each case, suggesting that there is only one structural isomer present. They can be assigned to increasing ring sizes, *i.e.* [12]-, [18]-, [24]- and [30]CPPs. The IMS trace of P_{36} , however, displays four clearly visible signals with significant drift time (DT) differences pointing towards, at least, four different isomeric structures. From the continuous increase in drift time with increasing CPP size of the homologous series in Figure 2, we assign the signal with the highest drift time to the macrocycle with 36 phenylene rings. The peaks at lower drift times, however, must be due to molecules with considerably more compact molecular shapes, *i.e.* smaller CCS values.




Figure 2. Ion mobility spectra of different CPP oligomers obtained at very low collision energy of 10V to remove non-covalently bonded aggregates (see SI).

Application of CID to the isomers of P_{36} after ion mobility separation provides insight into the individual fragmentation behavior (Figure S2). In Figure 3, the intact isomeric structures of P_{36} are denoted from M1 to M4 according to their drift time differences. The dissociation products own the same drift time as the precursor structures inside the green frames.




Figure 3. IMMS image of P_{36} obtained with CID at 200V after ion mobility separation. The x-axis corresponds to the m/z value and the y-axis represents the drift time. The brightness of the color indicates the peak intensity (yellow implies high intensity).

In case of M3, two dissociation products appear at 912 Da (annotated as P_{12}) and 1814 Da (annotated as P_{24}). These m/z values correspond exactly to [12]CPP and [24]CPP, therefore suggesting that M3 consists of two smaller macrocycles with different sizes. Analogously, there is only one intense peak at 1312 Da (P_{18}) with the same drift time as its parent ion M2. This indicates that M2 consists of two [18]CPPs. We thus safely conclude that the isomeric structures originate from catenanes with mechanically interlocked CPPs, where at least one covalent bond has to be broken to cause a dissociation into cyclic and linear fragment ions, both having exactly the same m/z values. M4 does not produce any fragments at 200 V which is in accordance with the assignment to the regular macrocycle [36]CPP. In case of M1, the smallest component of the P_{36}

isomer mixture, no dissociation into defined masses of smaller ring sizes can be observed. The absence of defined fragments and the compact structure (smaller than $2 \times [18]\text{CPP}$) suggest the presence of a molecular knot^[9] of the $[36]\text{CPP}$. Furthermore, it can be speculated that the high signal width in case of M1 and M3 arises from conformational flexibility of certain isomers as well as from contributions of more complicated interlocked structures of the catenanes.

In summary, the results from fragment ion analysis of P_{36} provide clear evidence for the existence of two different catenanes, exclusively composed of either one $[12]\text{CPP}$ and one $[24]\text{CPP}$ or two $[18]\text{CPP}$ s as presented in Figure 4. Estimated from the peak intensities of P_{36} isomers, the catenanes occur in significant amounts as compared to the regular $[36]\text{CPP}$.




Figure 4. Isomeric structures of macrocycle, catenanes and trefoil knot hidden in P_{36} ($\text{C}_{216}\text{H}_{144}$, exact mass = 2737.1323 Da).

The assignment to catenanes with different compositions can be rationalized by molecular modeling. Therefore, the three-dimensional (3D) structures of the four P_{36} isomers with different topologies are energy optimized by molecular dynamics calculations^[11] and submitted to CCS value calculations (Experimental part in SI). The derived CCS values provide a relative size determination which can be qualitatively compared to the measured drift times in IMS. Ions with larger CCS values at the same molecular weight move slower in ion mobility separation, hence having higher drift times. The calculated CCS values are summarized in Table 1.

The size decrease, represented by the decrease of CCS values in the order of $[36]\text{CPP}$ (904 \AA^2), $[2]\text{catenane}$ ($[12]\text{CPP} + [24]\text{CPP}$) (878 \AA^2), $[2]\text{catenane}$ ($2 \times [18]\text{CPP}$) (841 \AA^2) and the trefoil knot of $[36]\text{CPP}$ (655.01 \AA^2), is in line with the

compositions and the decreasing drift times detected from M4 to M1.

Table 1. Calculated CCS values for the (single chain) macrocycle and potential catenane as well as trefoil structures of P_{36} .

	$[36]\text{CPP}$	$[2]\text{catenane}$ ($[12]\text{CPP} + [24]\text{CPP}$)	$[2]\text{catenane}$ ($2 \times [18]\text{CPP}$)	Trefoil knot
CCS value (\AA^2)	904 (Figure 4a)	878 (Figure 4b)	841 (Figure 4c)	655.01 (Figure 4d)

Thus, ion mobility tandem mass spectrometry of an insoluble oligomeric mixture of CPPs reveals an unexpected, but exciting structural complexity. Whereas the results from traditional mass spectrometry would have been erroneously assigned to pure macrocycles up to very large ring sizes with 84 phenylene units, our findings demonstrate the co-existence of macrocycles and catenanes. Even though this has been proven for isomeric structures at a mass of 2737 Da, it can be expected that the larger macrocycles are in reality even more complex mixtures of catenanes and knots rather than ordinary monocycles.

From a synthetic viewpoint, the formation of significant amounts of CPP catenanes without targeted approaches of supramolecular preorganization^[12] of the phenylene building blocks is certainly surprising and deserves further attention.

Acknowledgements


This work was supported by the European Research Council grant on NANOGRAPH, Graphene Flagship (No. CNECT-ICT-604391), European Union Projects UPGRADE and MoQuaS (contract N.610449).

Keywords: Cycloparaphenylene • Ion mobility mass spectrometry • Catenane • Molecular Knot • Macrocycle

- [1] aR. Jasti, C. R. Bertozzi, *Chemical physics letters* **2010**, *494*, 1-7; bH. Omachi, T. Nakayama, E. Takahashi, Y. Segawa, K. Itami, *Nature chemistry* **2013**, *5*; cH. Omachi, Y. Segawa, K. Itami, *Accounts of chemical research* **2012**, *45*, 1378-1389.
- [2] aX. Guo, M. Baumgarten, K. Müllen, *Progress in Polymer Science* **2013**, *38*, 1832-1908; bC. Li, M. Liu, N. G. Pschirer, M. Baumgarten, K. Müllen, *Chemical reviews* **2010**, *110*, 6817-6855; cM. R. Talipov, A. Boddada, Q. K. Timerghazin, R. Rathore, *The journal of physical chemistry. C, Nanomaterials and interfaces* **2014**, *118*, 21400-21408.

- [3] E. R. Darzi, R. Jasti, *Chemical Society Reviews* **2015**, *44*, 6401-6410.
- [4] R. Jasti, J. Bhattacharjee, J. B. Neaton, C. R. Bertozzi, *Journal of the American Chemical Society* **2008**, *130*, 17646-17647.
- [5] aE. R. Darzi, T. J. Sisto, R. Jasti, *The Journal of Organic Chemistry* **2012**, *77*, 66246628; bF. E. Golling, M. Quernheim, M. Wagner, T. Nishiuchi, K. Müllen, *Angewandte Chemie International Edition* **2014**, *53*, 1525-1528; cM. Quernheim, F. E. Golling, W. Zhang, M. Wagner, H.-J. Räder, T. Nishiuchi, K. Müllen, *Angewandte Chemie International Edition* **2015**, *54*, 10341-10346; dF. E. Golling, S. Osella, M. Quernheim, M. Wagner, D. Beljonne, K. Müllen, *Chemical Science* **2015**, *6*, 7072-7078.
- [6] aL. Przybilla, J.-D. Brand, K. Yoshimura, H. J. Räder, K. Müllen, *Analytical chemistry* **2000**, *72*, 4591-4597; bK. Yoshimura, L. Przybilla, S. Ito, J. D. Brand, M. Wehmeir, H. J. Räder, K. Müllen, *Macromolecular Chemistry and Physics* **2001**, *202*, 215-222; cS. Trimpin, S. Keune, H. J. Räder, K. Müllen, *Journal of the American Society for Mass Spectrometry* **2006**, *17*, 661-671.
- [7] aA. B. Kanu, P. Dwivedi, M. Tam, L. Matz, H. H. Hill Jr, *Journal of Mass Spectrometry* **2008**, *43*, 1-22; bR. Mukhopadhyay, *Analytical chemistry* **2008**, *80*.
- [8] aG. De Bo, J. De Winter, P. Gerbaux, C.-A. Fustin, *Angewandte Chemie International Edition* **2011**, *50*, 9093-9096; bC. a. Scarff, J. R. Snelling, M. M. Knust, C. L. Wilkins, J. H. Scrivens, C. A. Scar, *Journal of the American Chemical Society* **2012**, *134*, 9193-9198.
- [9] C. A. Schalley, J. Hoernschemeyer, X. Li, G. Silva, P. Weis, *International Journal of Mass Spectrometry* **2003**, *228*, 373-388.
- [10] aA. Yol, C. Wesdemiotis, *Reactive and Functional Polymers* **2014**, *80*; bK. Guo, Z. Guo, J. M. Ludlow, T. Xie, S. Liao, G. R. Newkome, C. Wesdemiotis, *Macromolecular rapid communications* **2015**, *36*, 1539-1552; cW. Zhang, M. Quernheim, H. J. Räder, K. Müllen, *Analytical Chemistry* **2016**, *88*, 952-959.
- [11] M. D. Hanwell, D. E. Curtis, D. C. Lonie, T. Vandermeersch, E. Zurek, G. R. Hutchison, *Journal of cheminformatics* **2012**, *4*, 17-17.
- [12] aM. Weck, B. Mohr, J.-P. Sauvage, R. H. Grubbs, *The Journal of Organic Chemistry* **1999**, *64*, 5463-5471; bG. Gil-Ramírez, D. a. Leigh, A. J. Stephens, *Angewandte Chemie International Edition* **2015**, *54*, 6110-6150.

Solvent-free MALDI ion mobility mass spectrometry is used to analyze an insoluble product mixture of cycloparaphenylenes. For the first time, monocycles and isomeric catenanes can be detected in a simple reaction mixture, where only monocycles could be expected.



Wen Zhang, Ali Abdulkarim, Florian E. Golling, Hans Joachim Räder*, and Klaus Müllen*

Page No. 1– Page No.4

Cycloparaphenylenes and their Catenanes – Complex Macrocycles Unveiled by Ion Mobility Mass Spectrometry

WILEY-VCH