





# Synthesis of transition metal molybdates and tungstates and their reactivity in alkane and alkene oxidation

M. Jastak, R. Schlögl, A. Trunschke Fritz-Haber-Insitut der Max-Planck-Gesellschaft, Berlin, Germany

# Introduction


Ni-Co molybdates have been reported to show high selectivities for propene in the oxidative dehydrogenation of propane [1, 2]. Other binary transition metal molybdates of the type  $MTeMoO_{\theta}$  (M = Co) were reported to be active and selective in the partial oxidation of propene to acrolein and acrylic acid [3]. With the addition of Mn, sodium tungstate, which is structurally related to the corresponding molybdate, also showed considerable activity in the oxidative coupling of methane [4]. In this work, we have synthesized phase-pure  $MMoO_{\theta}$  and  $MTeMoO_{\theta}$  (M = Co, Mn) catalysts via a hydrothermal route and analysed their properties.

### Experimental

Phase-pure  $\alpha$ -MnMoO<sub>4</sub> was prepared by hydrothermal synthesis at 190°C for 20 h (13 bar pressure). Initially two solutions were prepared containing Mn(NO<sub>3</sub>)<sub>2</sub> x 4H<sub>2</sub>O and Na<sub>2</sub>MoO<sub>4</sub> x 2H<sub>2</sub>O. The latter solution had a Molybdenum concentration of 0.2 mol/l. Both solutions were mixed together at room temperature in the autoclave-vessel to form a white slurry (pH = 5.4). The suspension obtained after hydrothermal synthesis was filtered and washed with distilled water (150 ml). The solid fraction was dried at 100°C for 16 h.

# Results

The  $\alpha$ -MnMoO<sub>4</sub> phase prepared by hydrothermal synthesis was analyzed by X-ray diffraction (XRD), X-ray fluorescence (XRF), nitrogen adsorption, thermal analysis, temperature-programmed reduction (TPR) and electron microscopy. Figure 1 shows the X-ray diffraction pattern of the prepared  $\alpha$ -MnMoO<sub>4</sub> sample after drying. It is clearly observable that the  $\alpha$ -MnMoO<sub>4</sub>- phase is the only crystalline phase that is obtained after hydrothermal synthesis. Also no traces of Na have been found by XRF- analysis of this sample. The as-prepared  $\alpha$ -MnMoO<sub>4</sub>- phase shows a high thermal stability for temperatures of 500 °C and a quite high specific surface area of around 87 m²/g. The TPR- profile of the synthesized  $\alpha$ -MnMoO<sub>4</sub>-phase (see Figure 2)






Figure 1: X-Ray powder diffractogram of phase-pure α-MnMoO<sub>4</sub> prepared by hydrothermal synthesis

Figure 2: TPR-measurement of the synthesised  $\alpha$ -MnMoO<sub>4</sub>-phase (5.24% H<sub>2</sub> in Ar, 80ml/min)

shows an onset of H<sub>2</sub> consumption at about 410°C and two major peaks of H<sub>2</sub>-consumption with their maximums at about 590 and 880°C respectively. The onset and temperatures of maximum H<sub>2</sub> consumption are substantially lower than those reported for catalysts synthesized by other preparation routes [5]. This could be due to a higher oxygen mobility of the samples synthesized by the relatively fast and mild hydrothermal route.

### Conclusions

The here presented results show an example of a simple synthesis route to molybdate- and tungstate phases, which could potentially be used in selective oxidation reactions. The high thermal stability, the relatively high specific surface area, and increased oxygen mobility concluded from TPR-measurements of the synthesized sample suggest the applicability of hydrothermally synthesized molybdates and tungstates as catalysts in selective oxidation reactions of alkanes and alkenes.

#### References

- Y. -S. Yoon, N. Fujikawa, W. Ueda, Y. Moro-oka, K. -W. Lee; Catalysis Today 24 (1995), 327-333
- D. L. Stern, R. K. Grasselli; Journal of Catalysis 167 (1997), 550-559
- [3] B. Grzybowska, A. Mazurkiewicz, J. Sloczynski; Applied Catalysis 13 (1985), 223-243
- [4] Z. Gholipour, A. Malekzadeh, R. Hatami, Y. Mortazavi, A. Khodadadi; Journal of Natural Gas Chemistry 19 (2010), 35-42
- [5] S. A. Driscoll, L. Zhang, U. S. Ozkan; in "Catalytic Selective Oxidation" 523 (1993) 340-353 American Chemical Society, Washington DC